

# SomaLift

## The **Powered** Toilet Raiser








**Safety • Independence • Dignity**

The SomaLift provides powered assistance  
for standing up from the toilet

### Key Features:

- |                          |                        |
|--------------------------|------------------------|
| • powered lift mechanism | • easy to fit and move |
| • simple hand control    | • easy to use          |
| • use with a toilet      | • easy to clean        |
| • use as a commode       | • BUY or RENT          |

	<ul style="list-style-type: none"> <li>• Compact, lightweight design for easy fitting</li> <li>• Smooth, wipeable surfaces make cleaning easy</li> <li>• Non-slip surfaces and armrests enhance safety</li> </ul>
	<ul style="list-style-type: none"> <li>• Powered lift mechanism raises and lowers the seat</li> <li>• Simple controls in the armrest allow users to adjust the seat height effortlessly</li> </ul>
	<ul style="list-style-type: none"> <li>• A removable rechargeable <b>battery</b> is in the frame</li> <li>• Power indicator will tell you when to recharge</li> </ul>
	<ul style="list-style-type: none"> <li>• Removable <b>backrest</b> - offers comfort and security</li> <li>• Removable <b>commode</b> pan</li> <li>• Can be used as a standard toilet seat</li> </ul>
	<ul style="list-style-type: none"> <li>• <b>Height-adjustable</b> to fit various toilet heights</li> <li>• <b>4 legs</b> to allow use in a static position OR</li> <li>• <b>4 castors</b> (2 with brakes) for easy movement</li> </ul>

### 3 product packages available:

<b>GOLD</b>	<b>SILVER</b>	<b>BRONZE</b>
<ul style="list-style-type: none"> <li>• 2 batteries</li> <li>• Pre-assembled</li> </ul>	<ul style="list-style-type: none"> <li>• 2 batteries</li> <li>• Assembly required</li> </ul>	<ul style="list-style-type: none"> <li>• 1 battery</li> <li>• Assembly required</li> </ul>

4 Northgate, Crown Road, Enfield, London. EN1 1TG  
020 8805 2020 | [info@fortunabambach.com](mailto:info@fortunabambach.com)  
[www.fortunamobility.com](http://www.fortunamobility.com)

