

# WHEELCHAIRS

## Wheelchair Range:



**Transport** - 4 equal size wheels and are suitable for occasional appointments.



**Attendant (Transit)** - 4 small wheels and locking brakes on the handles. Our range is made from aluminium with a variety of configurations to suit the user's needs.



**Self-propelled** - larger wheels at the back with a rim that the user can grip to push the wheelchair forward. Made from aluminium with a variety of configurations to suit the user's needs.



## Specialist Wheelchairs:



**Tilt in Space** - tilt backwards whilst raising the legs for users to assist pressure relief and safety. Available in variety of models and sizes, from lightweight folding options which are ideal for transporting by car, or heavier wheelchairs which are set up for enhanced comfort and all day use.

**Bespoke (or made to measure)** - mainly used by those who are in their wheelchair for most of the day and are set up to the user's exact needs. The range of options is vast and due to their highly bespoke nature a thorough assessment by an expert is essential.



**Paediatric Wheelchairs** – we stock self-propelled and tilt in space wheelchairs in paediatric sizes but if something more specific is needed we can provide advice and arrange an assessment and have the wheelchair made to measure.



# WHEELCHAIRS

## Powerchairs

Powered wheelchairs can provide independence to those who are unable to use a self-propelled wheelchair.

Powerchairs can be bought in several forms to suit the person's individual needs:

- Indoor powerchairs are smaller and exceptionally manoeuvrable, to make it easier to move around in the house
- Outdoor powerchairs, are generally bigger, with larger wheels to maximise stability and comfort
- Bariatric powerchairs for larger users
- Folding powerchairs make it easier to transport by car



## Wheelchair Accessories

Our range of accessories improve both user and carer experience of the wheelchair.



To learn more visit our website:  
**[www.fortunamobility.com/wheelchair-overview](http://www.fortunamobility.com/wheelchair-overview)**

**FORTUNA.BAMBACH**  
Advice Choice Service

**[info@fortunabambach.com](mailto:info@fortunabambach.com) | [www.fortunamobility.com](http://www.fortunamobility.com)**

4 Northgate Business Centre, Crown Road, Enfield, London. EN1 1TG

020 8805 2020

