

ATTENDANT PROPELLED SHOWER COMMODE CHAIR

Our attendant Propelled Shower Commode Chair is manufactured from the highest quality materials with a rust-free alloy frame. A seamless PU seat and back provides maximum comfort and safety. The Shower Commode Chair features padded hinged arm supports for side transfers, a detachable padded back assembly and swing away footrests.



Specification

Product Dimensions (mm): 980x550x520
Net weight (kg): 14
Maximum User Weight (kg): 127
Maximum User Weight (stones): 20
Overall Depth (mm): 820 incl. footrests
Overall Height (mm): 910
Seat Depth (mm): 425
Seat Height (mm): 520
Seat Width (mm): 435
User Weight Limit (kg): 127

To learn more call us on
020 8805 2020 or visit our website :
www.fortunamobility.com/commodes
info@fortunabambach.com | www.fortunamobility.com
020 8805 2020

4 Northgate Business Centre, Crown Road, Enfield, London. EN1 1TG



FORTUNA.BAMBACH
Advice Choice Service

