

Seahorse Plus



3 in 1 Shower, Toileting & Commode Chair

Developed from the Seahorse Sanichair® it is designed to be supportive and safe. In addition to being highly adjustable it makes showering easier for both user and carer.



Seahorse Plus 3 in 1 Shower, Toileting & Commode Chair

The chair is made from high quality plastic and aluminium, so it is light, easy to clean and resistant to corrosion.

Available in 3 chair sizes, it is suitable for children from age 5 through to adults up to 1.7m, with the added potential to make further adjustments with each chair size.

The Seahorse Plus has a wide range of accessories (see overleaf) making it one of the most competitive showering & toileting chairs in the market.

Features

- Gas assisted tilt in space up to 40°
- Simple push button locking tilt for security & ease of use
- 100mm height adjustments
- Height adjustable flip-up footrests with foot-straps
- Height & depth adjustable headrest
- Lap belt
- Removable armrests
- Removable pommel
- Large castors with brakes for ease of movement and maximum security
- Designed for use over most standard toilets or with the removable potty

**To learn more or to arrange a joint visit contact us
on 020 8805 2020 or visit our website:**

www.fortunamobility.com/sea-horse



Chair Size Options

Small: age range guidance 5-10yrs, safe working load 64kg/10 stone

Medium: age range guidance 10-15yrs, safe working load 102kg/16 stone

Large: age range guidance 15-18yrs + smaller adults, safe working load 152kg/24 stone, maximum user height approx. 1.7m/5' 7"



Accessories

Liners

Removable, soft foam liners with interchangeable widths make the Seahorse Plus comfortable, adaptable and hygienic.

The liners have two separate back and seat sections for easy cleaning and for modular, bespoke fitting.

A choice of thickness permits greater precision in fitting the seat to the user and enables the chair to adapt to the growth of a child.



Anterior Supports

An option to help with positioning and to aid the toileting process.

The flip-up design allows easy entry and they can also be removed when not required.



Harness Options Chest Harness and Butterfly Harness options available for users needing more secure positioning.

**To learn more or to arrange a joint visit contact us
on 020 8805 2020 or visit our website:**

**www.fortunamobility.com/sea-horse
info@fortunabambach.com | www.fortunamobility.com**

4 Northgate Business Centre, Crown Road, Enfield, London. EN1 1TG



FORTUNA.BAMBACH
Advice Choice Service

