

**FITS  
PROFILING BEDS**

**EASY  
INSTALLATION**

**GREAT  
VALUE**

**FORTUNA.BAMBACH**  
Advice Choice Service



# EASYFIT PLUS+ Bed Rail



The EasyFit Plus+ Bed Rail provides a stable support for users getting into and out of bed and can also be used to assist the user when repositioning in bed.

It is designed for use with most types of beds including; divans with feet / castors, framed slatted beds and non-height adjustable profiling beds.

The sleek, sophisticated design of the EasyFit Plus+ Bed Rail makes it a complementary addition to any bedroom environment whilst promoting freedom and independence.

## FEATURES:

- **Fit most beds** including non-height adjustable profiling beds.
- **Quick and easy to install** – no tools or straps required.
- **Secures easily** - unique clamps slide quickly in to position under the bed.
- **Stable** - wide, low base plate and splayed legs.
- **Stores flat** when not in use.
- **Weight capacity** - 25st./160 kg
- **Standard Height model**



**To learn more or to arrange a joint visit contact us on  
020 8805 2020 or visit our website  
[www.fortunamobility.com](http://www.fortunamobility.com)**

[info@fortunabambach.com](mailto:info@fortunabambach.com) | [www.fortunamobility.com](http://www.fortunamobility.com)

020 8805 2020

4 Northgate Business Centre, Crown Road, Enfield, London. EN1 1TG



# Quick and Easy to Install

**1.** The EasyFit Plus+ attaches to its base without the need for tools. Simply slot the handrail over the base section and position EasyFit Plus+ under the bed.



**2.** Slide the unique clamps into contact with the underside of the bed frame – wide pads protect the underside of the bed.



**3.** Lock the clamps at the required height by closing the levers.



**4.** Turn the adjustment lever a couple of turns to firmly secure the rail to the bed frame.



## SPECIFICATION:

**Overall Width:** 522mm (20 ½")  
**Overall Length:** 712mm (28")  
**Overall Height:** 820mm (32 ¼")  
**Top Bar Height:** 692mm (27 ¼")

**Base Width:** 300mm (11 ¾")  
**Base Depth:** 180mm (7")  
**Handrail Width:** 285mm (11 ¼")

**FORTUNA.BAMBACH**  
Advice Choice Service

info@fortunabambach.com | www.fortunamobility.com  
020 8805 2020

4 Northgate Business Centre, Crown Road, Enfield, London. EN1 1TG