

StairSteady

For Confidence & Safety

The StairSteady helps you to keep your independence while staying safe. It's ideal for anyone who uses a walking aid or rollator:

Walking Stick



Stair Rail



Walker




StairSteady

Wheelchair



Stairlift

- ✓ Cost-effective: it can be less than HALF the price of a stairlift
- ✓ Keep mobile and active
- ✓ Feel safer and more confident
- ✓ Discrete and tidy within the home



FORTUNA.BAMBACH
Advice Choice Service



Call us to learn more

The StairSteady includes a:

- ✓ Rail which fixes to the wall
- ✓ Sliding handle that locks when weight is applied while going up or down stairs

With no motorised parts, the StairSteady is simple to use and extremely reliable.



Testimonial:

"I live in an upstairs flat and I am really pleased with my StairSteady. I feel much safer on the stairs because it feels so sturdy. I do go up and down more often now. The StairSteady is worth every penny and is much cheaper than a chair lift. Really the best thing I have purchased in years."

Mrs Tattershall



020 8805 2020 | info@fortunabambach.com

Book your **FREE** assessment with one of our technicians

www.fortunamobility.com/stairsteady

4 Northgate Business Centre, Crown Road,
Enfield, London EN1 1TG

