

# LIFT SEAT

**Powered toilet raiser** helps the user retain their independence and dignity in the bathroom. At the same time the toilet lift will help reduce effort and manual handling risks for carers who provide toileting assistance.



- Lowers and raises a user for using the toilet safely
- Utilises an anatomically correct "sit-to-stand" motion
- Use over your existing toilet, no conversion work required
- Oval-shaped opening for easy toileting
- Reduces the need for assistance
- Compatible with Closomat/ Gerberit

## Features:

1. **Handheld Control:** 2 button handheld control lowers or raises the user.
2. **Flat Transfer Seat:** makes transfer easy, and is smooth with rounded edges.
3. **Bedside Commode:** as the Lift Seat does not need ground fixing over the toilet nor anywhere else, it is quick and easy to transport and place in the required position
4. **Power Options:** to ensure continuity of use comes complete with 2 batteries, with one being plugged into the motor and a second on charge. They can then be swapped over to ensure continuity of use. A cable is also provided for hard wiring into a fuse spur, or in to a shower toilet, such as Geberit or Clos-O-Mat.
5. **Levelling Feet:** the 4 feet act as self-levellers if an uneven floor is encountered. The Lift Seat also has a range of 4 heights of feet that allow clients to adjust the height of the unit to better suit their needs.
6. **Streamlined Design:** the streamlined and compact design allows the Lift Seat to fit in virtually every type of bathroom and over every type of toilet.

**To learn more call us on**

**020 8805 2020 or visit our website :**

**[www.fortunamobility.com/bathroom](http://www.fortunamobility.com/bathroom)**

**[info@fortunabambach.com](mailto:info@fortunabambach.com) | [www.fortunamobility.com](http://www.fortunamobility.com)**

**020 8805 2020**

4 Northgate Business Centre, Crown Road, Enfield, London. EN1 1TG



**FORTUNA.BAMBACH**  
Advice Choice Service

