

# e-FOLDI

## Easy fold scooter



- Ideal for visiting shops or day trips
- Appropriate for both indoor and outdoor use
- Adjustable tiller for comfortable grip
- Pneumatic tyres
- 14 miles on single charge
- 12-month manufacturer warranty

### Technical Specification:

- Scooter length: 117 cm (46")
- Scooter width: 40.6 cm (16")
- Motor: 250 watts
- Maximum speed: 4 mph
- Battery size: 10 Ah
- Distance: 16km/ 10 miles \*
- Total weight: 15 kg
- Weight capacity: 120 kgs (19 stone)

\* Passenger weight, ground conditions, temperature and battery condition/size can affect travel distance



To learn more visit our website:  
[www.fortunamobility.com/wheelchair-overview](http://www.fortunamobility.com/wheelchair-overview)

**FORTUNA.BAMBACH**  
Advice Choice Service

[info@fortunabambach.com](mailto:info@fortunabambach.com) | [www.fortunamobility.com](http://www.fortunamobility.com)

4 Northgate Business Centre, Crown Road, Enfield, London. EN1 1TG

020 8805 2020

