



SomaCare ACER

A Modular Clinical High Back Chair

There has been limited choice of high back chairs,

- Seats with metal frames, which although inexpensive have non-breathable vinyl coverings and non-clinical foam.
- 'Fireside chairs' are attractive, as furniture – but do not address most clinical needs.
- Some healthcare chairs exist, but are non-adjustable and expensive.


Prioritising Healthcare Needs

Our Solution: The SomaCare ACER Range

Following extensive discussion with professionals, we've designed a range of seats at different price points, which can be selected for the healthcare needs of the user as they evolve.

Features

- All seat cushions have **breathable fabric** and **pressure-relieving foam**. According to the version chosen, this can also be specified for the backrest and armrests.
- A high back chair can make it easier to stand from sitting, while promoting the retention of function.
- Enhanced **adjustability** for seat height and depth as well as the configuration of the backrest.
- The Acer's **modular** design
 - increases choice when specifying initial requirements.
 - allows the specification to evolve with changing needs.
 - The UK-made seat is easy to repair and can be re-used for different clients, reducing the **carbon footprint**.
- The comfortable design will promote psychological well-being and calm
- The correct seating will reduce carer needs
- The design provides a highly **cost-effective** solution



Frame base for modular design.

AC1 This is the baseline model. The seat cushion is covered with breathable fabric and has medium-risk pressure foam. **£265 ex VAT**

FB-AC1-CSH

AC2 This seat has enhanced protection for the main pressure areas – the seat cushion and arms. **£295 ex VAT**

The seat cushion is covered with breathable fabric and has medium-risk pressure foam while the armrests are covered in Dartex and a layer of pressure-relieving foam.

FB-AC2-ARR

AC3 All contact areas are covered in breathable fabric and pressure-relieving foam. **£375 ex VAT**

The seat cushion is covered with breathable fabric and has medium-risk pressure foam while the armrests and backrest are covered in Dartex and a layer of pressure-relieving Vasco foam.

FB-A3-BCK

AC4 This seat adds a Dartex, shaped cushion (Vasco foam) **£395 ex VAT**

As well as providing an attractive design, this seat cushion can be reversed to provide greater depth adjustability. There is an option order arm infills for a feeling of greater security.

FB-A4-DTX

AC5 If a user is at high risk of pressure problems, this seat comes with a platform for use with a high-risk pressure cushion, including air-cushions, which can be specified by their healthcare team. **£330 ex VAT**

FB-A5-BAS



AC2



AC4



AC4 & Armrest Infills

Choice of Colours

Discounts available for larger orders

www.fortunamobility.com/blog

support@fortunabambach.com or **020 8805 2020**

Fortuna Mobility Centre, 4 Northgate, Crown Road, Enfield, EN1 1TG

