

A photograph showing a person lying in a bed, covered with a dark blue blanket. The bed is equipped with the Hugga Slide system, which consists of several large, dark brown, rectangular cushions or pads positioned around the person's body to facilitate movement and reduce pressure. The person is wearing blue trousers and a dark top. The background is a plain, light-colored wall.

How do we cope with Pressure
injuries, Single handed care & Posture
in **one piece of equipment?**

with the new **Hugga Slide**

In collaboration with **Careability Healthcare & Parafricta** low friction technology, **Posture Care** has developed the "**Hugga Slide**" 4 way slide, in bed management system.

- Designed to help carers perform moving and handling tasks as well as postural management in bed
- The base sheet has a low friction panel which can cover the entire bed length
- High friction side panels to prevent the user from sliding down the bed
- Top sheet is compatible with the Hugga sleep system www.hugga.co.uk for night time positioning in bed
- Reposition the user in all directions
- Achieve midline without heavy lifting
- Safer, easier turns from supine to side lying
- Sheets can be washed up to 80 degrees for infection control

Contact us at enquiries@hugga.co.uk for more information
or a **free** demonstration/assessment